

Essential Basic

Quality Foam Mattress




A cost effective solution for comfort
and pressure redistribution






The new **Essential Basic** mattress offers the user a comfortable and supportive surface. This is an ideal entry level mattress for individuals who do not require therapy for pressure ulcer prevention.

- ▶ Multi-stretch polyurethane fabric cover
- ▶ Zipped around 3 sides for ease of inspection and laundering



Technical data ▼

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	 Size range ▼	 Max. user weight ▼	 Max. cleaning temperature ▼	Nominal density range (kg/m ³) ▼	Article no. ▼
Essential Basic	2000 x 880 x 150 mm	140 kg	95°	32 - 34	1586660

Mattresses is manufactured to:
EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

1. *Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement).*

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland