



Invacare® **Bed Accessories**

To fit all needs!



The *Invacare* bed range features a wide range of accessories which gives you the possibility to customise your bed both in terms of design and taste. This allows perfect fittings to personal environment with supporting accessories which tailors the bed to individual client needs.

In the catalogue you will also find compatibility references for beds, side rails and bed ends.

Bed ends



Victoria bed end

Detachable bed end for **SB 755**, **Alegio_{NG}** and **SB 400** in beech, white beech, cherry and light oak. Compatible with Line and Britt V full length side rails.



Sophie bed end

Available in beech, white beech, cherry and light oak. Standard/low version for **SB 755**, **Alegio_{NG}** and **SB 400**, compatible with Line and Britt V full length side rails.



Susanne bed end

Available in beech, white beech, cherry and light oak for **SB 755**, **Alegio_{NG}** and **SB 400**, compatible with Line and Britt V full length side rails.



Anita bed end

Available in beech, white beech, cherry and light oak for **SB 755**, **Alegio_{NG}** and **SB 400**, compatible with Line and Britt V full length side rails.



Emma bed end

Detachable bed end for **SB 755**, **Alegio_{NG}** and **SB 400** in beech, white beech, cherry and light oak. Compatible with steel side rails.



Piggy bed end

Detachable bed end in beech, white beech, cherry and light oak. Available for **SB 755**, **SB 400** and **Alegio_{NG}** and compatible with steel side rails.



Olivia bed end

Detachable bed end in beech and light oak, compatible with wooden side panels and Scala Basic Plus 2 side rail. Available for **SB 755**, **SB 400** and **Alegio_{NG}**.



Sabrine bed end

Available in cherry and beech for **Sonata** and **Alegio_{NG}**. Compatible with steel side rails.



Vibeke_{NG} bed end

Available in cherry for **SB 755**, **SB 400**, **Alegio_{NG}** and **Sonata**. Compatible with steel side rails.



Nadine bed end

Bed end with curved metal tubing and wooden panel in beech for **Sonata**. Compatible with steel side rails.



Valeria bed end

Bed end with curved metal tubing and wooden panel in beech for **Sonata**. Compatible with steel side rails.



Camila bed end

Bed end in cherry for **Sonata** and **Alegio_{NG}**. Compatible with the Flora and Lisa side rails.

Side panels



Wooden side panel

Available in beech, white beech, cherry and light oak.



Wooden side panel

Can be used in combination with Scala Basic Plus 2 side rail. Available in beech and light oak.

Side rails



Torill Long

Split side rails offer full protection when placed in their highest position. Available in beech, white beech and light oak.



Torill Short

Split side rails with opening in the middle. Available in beech, white beech and light oak.



Line

Full length aluminium side rail. Available in beech, white beech, cherry and light oak.



Line Extendable

Full length aluminium side rail. Available in beech, white beech, cherry and light oak. Adjustable length 200 - 220 cm.



Britt V

Full length wooden side rail in lengths of 200 cm, 210 cm and 220 cm. Available in beech, cherry, white beech and light oak.



Side rail height extender

Height of full length side rails can be increased by 15 cm.



Bella wooden side rail

Full length side rail in wood. Available in beech.

Side rails



Aria steel side rail

Full length side rail in steel.



Verso II side rail

A collapsible steel side rail.



Flora and Lisa wooden side rails

Full length wooden side rail available in cherry.



Diana siderail

A lowerable metal side rail.



Invacare *Scala Basic 2*

Collapsible steel side rail.



Invacare *Scala Basic Plus 2*

Collapsible steel side rail with plastic inserts at the mounting brackets.



Invacare *Scala Medium 2*

Collapsible steel side rail for higher mattresses.



Invacare *Scala Decubi 2*

Collapsible steel side rail for extra high mattresses.

Side rail covers



Cover for Line and Britt V side rails

Padded covers for Line and Britt V side rails. Cover remains in place when side rail is lowered. Available in padded, extra padded and extended + 10 cm and + 20 cm.



Net cover for Line and Britt V side rails

Net covers for Line and Britt V side rails. Cover remains in place when side rail is lowered. Available as standard and extended + 10 cm and + 20 cm.



Cover for metal side rail

Padded cover for metal side rails. Cover remains in place when side rail is lowered. Available in standard, extended, and extra padded.



Cover for side rail height extender

Padded cover for side rail height extender placed on Britt V/Line.

Support handles



Support handle

A very stable support handle for getting in and out of bed:

H25 X W30 cm / W80 cm

H40 X W30 cm / W50 cm / W95 cm



Swivel support handle

The swivel support handle makes it possible to swing the legs out over the handle before the grip is folded out. Can be used alone or in pairs.



Support handle for Octave

Makes getting in and out of bed much easier for users.

H40 x W30 cm.



Sit-up support

Help for getting in and out of bed.

For **SB 755**, **SB 400** and **Etude Plus**.



Stand-up support

Makes getting in and out of bed easier. Easy to use and mount onto the bed. Fits all Invacare beds.

Lifting poles



Lifting pole

Catching handle can be adjusted in both height and depth.



Swivel lifting pole

The lifting pole swings away to the side of the bed. Catching handle can be adjusted in both height and depth.



Lifting pole

Catching handle can be adjusted in both height and depth. For use with **Octave** beds.



Free-standing lifting pole


Height adjustable, catching handle can also be adjusted.

Other accessories




Distance bumper

Prevents walls from being damaged when the bed is moved. Can be placed on the base frame at the end of the bed or at the side of the bed.



Bumper wheels

The bumper wheel can be mounted horizontally or vertically. Available in three different sizes.



Castors

Different sizes of castors. Comes in 75, 80, 100, 125, and 150 mm.




Lowering kit

By using 75 mm wheels, **SB 755** can be lowered to 25 cm.



Cover tubes

To be mounted on the bed frame to cover the tubes when the bed is extended. For **SB® 755** and **SB® 400** in different sizes and colors.



Limitation crossbar

Prevents the leg support from dipping below horizontal position. Can be placed on the base frame. Available in ivory and aluminium. For **SB 755**.



Extension parts for Etude Plus

Available in 10 cm or 20 cm extension.



Extension for Octave

Extension for **Octave** for greater flexibility. 10 cm to 20 cm extension.



Mattress support extension for Alegio^{NG}

Extends the length of **Alegio^{NG}** by additional 15 cm.



Mattress support extension for Medley Ergo

Extends the length of **Medley Ergo** by additional 15 cm.



Wooden adapters - Octave bed

Wooden adapters for mounting **Octave** side rail, comes in beech and cherry.



Side rail lock

Strap for locking the side rail in its lowest position. Can be used with any side rail.

Other accessories



Transportation and storage rack

To be used with **SB 400** when the bed is dismantled.



Transportation fitting for *Etude Plus*

Fittings for easy transport of the **Etude Plus** beds.



Transportation rack

Can be used with **SB 755**. For transporting or storing an assembled bed.



Transportation rack

To be used for transporting a bed.




Transport trolley for *Alegio_{NG}*

Trolley for easy storage and transport. Can be dismantled.




Transport kit for *Alegio_{NG}*

Designed for use when storing and transporting the bed fully assembled.



Transport fitting for *Medley Ergo*

For easy storage and transport,



Bed Cradle

Bed cradle can be placed on the bed end to alleviate pressure from the duvet on the legs/feet.




IV Rack

Free-standing IV rack. Height-adjustable, castors for easy movement.



IV Rack

Can be placed in the lifting pole tube or on the bed frame.



IV Rack mount fittings

Fittings for mounting the IV rack either in the lifting pole tube or on the bed frame.

Electrical accessories



Hand control

Hand control for **Etude Plus** available without tilt, with anti-tilt only and both tilt and anti-tilt, with or without lockable functions.



Hand control

Hand control for **Alegio^{NG}** available with or without lockable functions.



Hand control

Hand control for **Medley Ergo** and **Accent**, available without tilt, anti-tilt only and both tilt and anti-tilt, with or without lockable functions.



Hand control

Hand control for **Sonata**.




Rheumatic hand control

Hand control for people with rheumatism



Holder for hand control

Flexible holder for hand control, easier access for users. Available for soft control, HB70, HB80, HL80 and HL85 and Dewert hand control. Comes in aluminium and ivory.



Out of bed alarm

This alarm fits most beds. Alarm can be adjusted to individual needs. Available in a standard and epileptic version.



LED lamp

Attractive, modern design, very low power consumption (4W). Can be mounted on the lifting pole bracket on all Invacare beds. Comes in aluminium.



IR Control

Infrared remote control for use when standard hand control is difficult to operate. Can be used for **SB 755** and **SB 400**.



IR Receiver (without transmitter)

The infrared receiver allows multiple functions to be controlled by the same remote control, e.g. bed, light and TV.



T-cable

Can be used when two hand controls are needed.

Electrical accessories



Battery back-up

To be used in case of power failure. The bed can be delivered prepared for battery back up. The battery back up will be installed afterwards.



Extension cable

2.5 metre extension cable for Linak hand control. Can be used for **SB 755** and **SB 400**, **Etude Plus** and **Octave**.



ACP box

Gives healthcare personnel the ability to lock the function of one or more motors. Can be used for **SB 755** and **SB 400**, **Etude Plus** and **Octave**.



Field power supply

Can be used to adjust the height or any other motor when no power is available. For service personnel only, not to be used with a client in bed.



Hand control

Hand control HB80 and HL80 with green button function. Each function of the HL80 can be individually locked. Available for **SB 755**.



Hand control

Hand control HB70. With green button function for user comfort. Available for **SB 400**.



Soft control

Function-based hand control with easy activated buttons. Available for beds with three and four motors.



Hand control for Octave

Hand control HL85. Each function on the HL85 can be individually locked.



Bed tray

Stable table/tray, fits directly on the side rail. Available in 85 and 90 cm.



Sonnet Basic

A stable working and reading table. Height-adjustable table with plate that can be adjusted from 65 - 100 cm in height. The work top can be angled in six different positions.



Pausa

Height-adjustable table with plate that can be adjusted from 72 - 116 cm in height. Max load: 10 kg.



Kauma

Height-adjustable table with plate that can be adjusted from 71 - 113.5 cm in height. Max load: 10 kg.



Sonnet Dual

Bedside table in beech, white beech, cherry and light oak with a practical top shelf, two drawers and small cupboard. The flap can be adjusted from 75 - 105 cm in height and set at an angle to suit user need.



Sonnet Combi

Bedside table in beech, white beech, cherry and light oak. All horizontal surfaces are in plastic laminate. Practical top shelf, one drawer and a cupboard. The flap can be adjusted from 31 - 93 cm in height.



Tendance

Bedside table in light beech, with one drawer and one cupboard.



Bedco Anti-Sliding

Prevents the backrest from sliding when the backrest is being operated.



Back Support

Backrest support can be used in an ordinary bed for added user comfort.



Back Support

Manual backrest support, with possible head support. Can be adjusted from a 45 - 60 degree angle.

Invacare **SB[®] 755, Alegio[™] NG, SB[®] 400**

Side rails	Bed ends				
	Victoria Sophie Susanne Anita	Sophie Low ² Susanne Low ² Anita Low ²	Emma Piggy VibekeNG ²	Camila (only for Alegio NG)	Olivia
Torill Long and Torill Short.	x	x			x
Line / Line Extendable / Britt V / Britt V Extended ¹	x	x			
Scala Basic 2 ¹ / Scala Basic Plus 2 ¹	x	x	x		x
Scala Medium 2 ¹	x	x	x		x
Scala Decubi 2 ¹	x	x	x		x
Verso II (w. bracket for 50x20 tube)	x	x	x		
Lisa				x	

Invacare **Sonata[®]**

Side rails	Bed ends	
	VibekeNG, Nadine, Sabrina, Valeria	Camila
Verso II	x	
Diana	x	
Flora	x	x

Invacare **Etude[®] Plus (Low), Medley[®] Ergo, Accent[™]**

Side rails	Etude Plus	Medley Ergo	Accent
Line / Line Extendable / Britt V / Britt V Extended ¹	x		
Diana			x
Scala Basic 2	x	x	x
Scala Basic Plus 2	x	x	x
Scala Medium 2	x	x	x
Scala Decubi 2	x	x	x
Full length side rail, Bella wood or Aria steel		x	x
Verso II (50x20)		x	x

Invacare **Octave[®]**

Side rails	Bed ends
Line	Fits Octave bed end
Line Extendable	
Britt V	
Britt V Extendable	

1. Versions do not fit with **Alegio_{NG}**
2. Only for 90 cm

Complies with IEC 60601-2-52

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International Sàrl
Benkenstrasse 260
4108 Witterswil
Switzerland

QR
Code

