



Invacare®

Essential Range



Invacare **Essential Basic**






The **Essential Basic** mattress is an ideal entry level mattress that offers a combination of comfort, pressure reduction and value for money. It features a high quality multi-stretch vapour-permeable, smooth and water resistant, polyurethane, fabric cover that helps to reduce shear and friction forces, in addition to helping meet modern hygiene standards.

- Multi-stretch vapour-permeable, water resistant, polyurethane, fabric cover
- Zip around three sides for ease of inspection and cleaning
- Fire retardant

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Size	Max. user weight	Max. cleaning temperature	Nominal density range (kg/m ³)	Article no.
					
Essential Basic	2000 x 880 x 150 mm	140 kg	80°	32 - 34	1586660

Invacare **Essential Care**






The **Essential Care** mattress is a pressure redistributing support surface, developed for home care and residential care environments. With the same high quality cover as the **Essential Basic**, the **Essential Care** mattress has a specially cut castellated foam surface, providing a cost effective solution for comfort and pressure redistribution to meet the needs of individuals 'At Risk/High Risk'¹ of developing pressure ulcers.

- Cut foam castellations
- Fire retardant



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Size	Max. user weight	Max. cleaning temperature	Nominal density range (kg/m ³)	Article no.
					
Essential Care (150 mm height)	2000 x 880 x 150 mm	140 kg	80°	32 - 34	1588619
Essential Care (120 mm height)	2000 x 880 x 120 mm	110 kg	80°	32 - 34	1588614

The **Essential Visco** mattress is a pressure redistributing support surface for the home care, residential and long-term care environment.

This mattress features a high quality multi-stretch vapour-permeable, smooth and water resistant, polyurethane, fabric cover that helps to reduce shear and friction forces on the delicate skin of the patient, in addition to helping meet modern hygiene standards.

The high density, visco elastic foam top layer conforms to the body shape offering maximum support and pressure relief.

This mattress meets the needs of individuals 'At Risk/High Risk' 1 of developing pressure ulcers.

- Visco elastic foam top layer
- Sturdy foam base for additional support
- Fire retardant



Features and options



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Size	Max. user weight	Max. cleaning temperature	Nominal density range (kg/m ³)	Article no.
Essential Visco (150 mm height)	2000 x 880 x 150 mm	140 kg	80°	32 - 34, Visco layer 48 - 52	1588533
Essential Visco (120 mm height)	2000 x 880 x 120 mm	110 kg	80°	32 - 34 Visco layer 48 - 52	1588538

The *Invacare Essential Plus* is an entry level pressure ulcer prevention mattress designed for the community healthcare environment. Practical and durable, this mattress offers comfort, support and pressure redistribution at a cost effective price. The **Essential Plus** is suitable for those considered at 'High Risk'¹ of developing pressure ulcers.

Effective Pressure Ulcer Prevention

The core of the *Invacare Essential Plus* is constructed from a series of 'keyhole cut' cells, which move independently to evenly distribute weight, promoting greater comfort and facilitating pressure ulcer prevention.

Cover is resistant to fluid ingress

The outer cover is a multi stretch, water resistant, polyurethane fabric providing excellent recovery and durability, which helps to reduce shear and friction forces. The cover also has a white inner liner allowing easy inspection for strikethrough or other cover damage.



Features and Options



Turn and rotate

Double sided castellated foam allows the mattress to be turned or rotated for increased product longevity.



Available in the following sizes:

L2000 x W855 x D140mm
Order code: 1576201

L1950 x W880 x D140mm
Order code: 1576200






Zip around three sides

Allows easy inspection of the foam core and easy removal for laundering.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Size	Max. user weight	Max. cleaning temperature	Nominal density range (kg/m ³)	Article no.
					
Essential Plus (195 mm length)	1950 x 880 x 140 mm	140 kg	80°	40 - 42	1576200
Essential Plus (2000 mm length)	2000 x 855 x 140 mm	140 kg	80°	40 - 42	1576201

1. Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement).

Mattresses is manufactured to: EN 597-1 Cig, BS EN 597-2 Match and BS 7177: Crib Source 5

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland

QR
Code

