

Invacare® | **Softform Excel**



Warranty
Foam 2 Years
Cover 2 Years



Pressure redistributing mattress

The *Invacare Softform Excel* provides a high level of comfort and is constructed from high specification foam which provides even weight distribution for greater comfort and pressure redistribution. With a robust foam base, it prevents the mattress from bottoming out and is suitable for those considered 'At Risk'* of developing pressure ulcers.



Yes, you can.®



Features and Options

Effective pressure ulcer prevention

The core of the *Invacare Softform Excel* is constructed from a series of 'keyhole cut' cells, which move independently to evenly distribute weight, promoting greater comfort and facilitating pressure redistribution. The robust base features a firm layer of high density foam to provide substantial support.

Multi stretch PU waterproof cover

The cover fabric provides excellent recovery and durability, and helps to reduce shear and friction forces. The base of the mattress is a tough coated polyurethane (PU) which protects the mattress when used with electric profiling bed frames.



Zip around three sides

Allows easy inspection of the foam core.



Rotate-only mattress

Rotate-only turning cycle will increase product longevity. Non turn reduces the amount of manual handling required.



Strong, durable construction

The *Invacare Softform Excel* is an established choice of healthcare professionals.

Technical data

Softform Excel	Standard width: 800, 855, 880, 900 mm	Standard length: 1950 mm	Height: 152 mm	Weight: 10 kg	Max. Client Weight: 140 kg (22 stone)

	Softform Excel	Cleaning: 80°
--	-----------------------	---------------

Manufactured to European Standards:
BS7177: Crib 5

Invacare Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2014 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Softform Excel - UK - 09/2014
1506260

Colour



* Essential nursing care is pivotal in pressure ulcer prevention. This mattress will positively contribute to the outcome of a pressure ulcer prevention care plan. Education, clinical judgement and action based planning based on vulnerability are fundamental factors in the prevention of pressure ulcers. A range of assessment scales can be used as a formal method of assessing risk from pressure ulcer development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

